



**PRODUCED IN
TRANSPARENT POLYMER**

SAME MATERIAL USED
IN AESTHETIC BRACKETS

- ✓ HIGHLY RESISTANT POLYMER
- ✓ HIGH RESISTANCE TO DISPLACEMENT
- ✓ SAFE AND RELIABLE REMOVAL
- ✓ HIGH MECHANICAL RESISTANCE TO BREAKAGE AND DISTORTIONS

Suggested Use:

Use light-curing resin (for example: Light Bond) for bonding the buttons;
Use liquid sealant at base of button;
We do not recommend the use of resin from single or chemically activated paste.



MAXIMUM TRANSPARENCY
AESTHETICS

TRAPEZOIDAL



OVAL



ROUND



Ideal combination with aligners

It is applied in malocclusion treatments dental, in movements such as gyroversions, tractions, anchoring systems and retracts. Ideal for aesthetic treatments.

Used with chain elastics, ligature wire and ligatures.

Available in formats: oval, round and trapezoidal.



1-877-889-6913 • (813) 263-3038

contact@aditekortho.com • www.aditekortho.com